
Test Machines Australia TMA-MC-1 Friction Coefficient Tester





Test Machines Australia



Main Application

The equipment is suitable for testing the dynamic and static friction coefficient of plastic film, sheet, paper and other materials. By measuring the smoothness of the material, we can control the opening of the packaging bag, the packaging speed of the packaging machine and other production quality process indicators to meet the requirements of the product.

Function and Characteristic

- 1. It meets GB, ISO, and ASTM test standard.
- 2. The flexible drive system drives the operation more smoothly and the test precision is higher.

Test Machines Australia

3. LCD full figure according to the test data, the results, curve

4. The quantitative and discrete analysis of the test results between the slippage of

the material and the test results are obtained.

5. Dynamic friction test, static friction test, free choice of static and static friction test,

flexible test.

6. The results of single piece, group test results were analyzed, and the output of

multiple reporting modes was micro-typed.

7. It can save 6 test data and results, with curve display, query and other necessary

functions.

Main Technical Parameters

load range: $0\sim5$ N

Accuracy: 0.5 grade

Travel: 70 mm 150 mm

Weight of the slider: 200g (Standard)

Test Speed: 150mm/min

Power requirements and dimensions: AC 220V, 50Hz, 400(L)×300(W)×180(H)mm

For further information please contact:



0418 369 505
sales@testmachines.com.au
www.testmachines.com.au

Test Machines Australia